

# HeartSine® samaritan® PAD

## family



	HeartSine samaritan PAD SAM 350P	HeartSine samaritan PAD SAM 360P	HeartSine samaritan PAD SAM 500P with CPR Advisor™
ICG-based real-time feedback on CPR rate and force ("Push harder," "Push faster," "Push slower," "Good compressions")			X
CPR coaching with metronome	X	X	X
Fully automatic shock delivery		X	
Semi-automatic shock delivery	X		X
High level of protection from dust and water ingress (IP56)	X	X	X
Integrated battery and electrode for one-change maintenance	X	X	X
4-year electrode and battery life	X	X	X
Low energy escalating protocol (SCOPE)	X	X	X
Operating impedance range of 20 – 230 ohms	X	X	X
Audible and visual metronome for CPR rate guidance	X	X	X
Compact, portable, and lightweight	X 1.1 kg (2.4 lb)	X 1.1 kg (2.4 lb)	X 1.1 kg (2.4 lb)
Pediatric compatible	X	X	X*
90-minute ECG recording time	X	X	X
Automatic self-test	X	X	X
Field upgradeable software	X	X	X
Carry case included	X	X	X

\* When used with Pediatric-Pak, CPR Advisor is disabled.

All claims valid as of 02/2025.


For further information, please contact your Stryker representative or visit our website at [stryker.com](https://www.stryker.com)

## Emergency Care Public Access

AED users should be trained in CPR and in the use of the AED. Although not everyone can be saved, studies show that early defibrillation can dramatically improve survival rates. AEDs are indicated for use on adults and children. AEDs may be used on children weighing less than 25 kg (55 lb) but some models require separate defibrillation electrodes.

The information presented is intended to demonstrate Stryker's product offerings. Refer to operating instructions for complete directions for use indications, contraindications, warnings, cautions, and potential adverse events, before using any of Stryker's products. Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your representative if you have questions about the availability of Stryker's products in your area. Specifications subject to change without notice.

Stryker or its affiliated entities own, use, or have applied for the following trademarks or service marks: CPR Advisor, HeartSine, Pad-Pak, Pediatric-Pak, samaritan, SCOPE, Stryker. All other trademarks are trademarks of their respective owners or holders.

 HeartSine samaritan PAD is a class III – 0123 device in accordance with EU MDR. Pad-Pak and Pediatric-Pak are CE class IIb – 0123 also in accordance with EU MDR.

 HeartSine samaritan PAD: UL Classified. See complete marking on product.

Date of issue: 02/2025

EN-UK

Made in U.K.

HeartSine SAM 500P is not available for sale in the U.S.

EC-HS-COMP-1535338\_REV-0\_en\_gb

Copyright © 2025 Stryker.



HeartSine Technologies, Ltd.  
207 Airport Road West  
Belfast  
Northern Ireland  
BT3 9ED  
United Kingdom  
Tel +44 28 9093 9400  
[heartssinesupport@stryker.com](mailto:heartssinesupport@stryker.com)  
[heartssine.com](https://www.heartssine.com)