

HeartSine samaritan® PAD 500P with CPR Advisor

Zoll AED Plus Product Comparison



Heartsine samaritan® PAD 500P Technical Overview



Physical	With Pad-Pak™ Inserted	
Size:	20 cm x 18.4 cm x 4.8 cm (8.0 in x 7.25 in x 1.9 in)	
Weight:	1.1 kg (2.4 lbs) including Pad-Pak Battery	

Energy Selection			
Adult:	Shock 1: 150J;	Shock 2: 150J;	Shock 3: 200J
Pediatric:	Shock 1: 50J;	Shock 2: 50J;	Shock 3: 50J

Defibrillator	
Waveform:	Self-Compensating Output Pulse Envelope (SCOPE) Biphasic waveform. Optimized biphasic escalating waveform compensates energy, slope and envelope for patient impedance

Charging Time	
New Battery:	Typically 150J in < 8 seconds, 200J in < 12 seconds
After 6 Discharges:	Typically 150J in < 8 seconds, 200J in < 12 seconds

Patient Analysis System	
Method:	Evaluates patient's ECG, signal quality, electrode contact integrity and patient impedance to determine if defibrillation is required
Sensitivity/Specificity:	Meets ISO 60601-2-4 and AAMI DF80

Event Documentation	
Type: Internal Memory	
Memory Capacity:	90 minutes of ECG (full disclosure) and event/incident recording
Playback Capabilities:	Custom USB cable directly connected to PC and Saver™ EVO Windows-based data review software

Environmental	
Operating/Standby Temperature:	o°C to 50°C (+32°F to +122°F)
Temporary Transportation Temperature:	–10°C to 50°C (14°F to 122°F) for up to two days. Unit must be returned to standby/operating temperature for 24 hours before use.
Relative Humidity:	5% to 95% (non-condensing)
Water Resistance:	IEC 60529/EN 60529 IP56
Altitude:	0 – 4,575 meters (0 to 15,000 feet)
Shock:	MIL STD 810F Method 516.5, Procedure I (40 G's)
Vibration:	MIL STD 810F Method 514.5+ Category 4 Truck Transportation – US Highways Category 7 Aircraft – Jet 737 & General Aviation (Exposure)
EMC:	EN 60601-1-2
Radiated Emissions:	EN55011
Electrostatic Discharge RF Immunity:	EN61000-4-3 80MHZ-2.5GHZ (10 V/m)
Magnetic Field Immunity:	EN61000-4-8 (3 A/m)
Aircraft:	RTCA/DO-160F, Section 21 (Category M)
Falling Height:	1 meter (3.3 feet)

Materials Used		
samaritan® PAD 500P:	ABS, Santoprene. Printed circuit board with electronic components.	
Housing:	ABS – Electrodes: Hydrogel, Silver, Aluminium and Polyester	

Pad-Pak — Electrode and Battery Cartridge Adult Pad-Pak (Pad-Pak-01) and Pediatric Pad-Pak (Pad-Pak-02)		
Shelf Life:	4 years from manufacture date	
Weight:	0.2 kg (0.44 lbs)	
Size:	10 cm x 13.3 cm x 2.4 cm (3.93 in x 5.24 in x .94 in)	
Battery Type:	Lithium Manganese Dioxide (LiMnO2)	
Capacity:	> 60 shocks at 200J 18V, 1.5 Amp Hrs	
Electrodes:	HeartSine samaritan® disposable defibrillation pads are supplied as standard with each device	
Placement:	Anterior-lateral (Adult); Anterior-posterior (Pediatric)	
Active Gel Area:	100 cm² (15.5 in²)	
Cable Length:	1 meter (3.3 feet)	

Lifesaving, Pure and Simple

U.S./Americas

HeartSine Technologies, Inc. 121 Friends Lane, Suite 400 Newtown, PA 18940 Toll Free: (866) 478 7463 Tel: +1 215 860 8100 Fax: +1 215 860 8192 info@heartsine.com

Europe/Rest of the World

HeartSine Technologies, Inc. 203 Airport Road West Belfast, Northern Ireland BT3 9ED Tel: +44 28 90 93 94 00 Fax: +44 28 90 93 94 01 info@heartsine.co.uk



The products described in this brochure all meet the applicable European Medical Directive requirements.

© 2014 HeartSine Technologies, Inc. All rights reserved. H009-014-124-0



practitioner.

