

HeartSine samaritan® PAD 500P with CPR Advisor

CardiaTech CardiAid® Product Comparison

Integrated
battery and
electrodes



HeartSine samaritan® PAD 500P

- ✓ **ICG-based CPR feedback**
- ✓ Compact, portable, lightweight (2.4 lbs/1.1 kg)
- ✓ Highest level of protection from dust and moisture ingress in the industry (IP56)
- ✓ **Integrated battery and electrodes**
- ✓ 4-year electrode life
- ✓ One item to track and replace after use or upon expiration
- ✓ **10-year warranty**
- ✓ Charging time typically 150J in < 8 seconds, 200J in < 12 seconds
- ✓ 90-minute ECG recording time
- ✓ Low escalating energy protocol (SCOPE™)
- ✓ Field upgradeable to latest AHA/ERC guidelines
- ✓ Pediatric compatible
- ✓ Audible and visual prompts for ease of use in a noisy environment
- ✓ Automatic self test
- ✓ Carry case included
- ✓ **Lower total cost of ownership; acquisition and operating expense**



CardiaTech CardiAid®

- ✗ No CPR feedback
- ✗ Larger, heavier (6.8 lbs/3.1 kg)
- ✗ Lower water resistance rating (IPX4)
- ✗ Separate battery and electrodes (Device must be returned to dealer every two years to replace battery)
- ✗ 3-year electrode life
- ✗ Two items to track and replace after use or upon expiration
- ✗ 6-year warranty
- ✗ Charging time < 24 seconds
- ✗ 240-minute ECG recording time (Recorded information can only be obtained by an authorized dealer)
- ✗ Fixed width energy protocol (Not FDA cleared)
- ✗ Not field upgradeable
- ✗ Not pediatric compatible
- ✓ Audible and visual prompts for ease of use in a noisy environment
- ✓ Automatic self test
- ✓ Carry case included

U.S./Americas

HeartSine Technologies, Inc.
121 Friends Lane, Suite 400
Newtown, PA 18940
Toll Free: (866) 478 7463
Tel: +1 215 860 8100
Fax: +1 215 860 8192
info@heartsine.com

Europe/Rest of the World

HeartSine Technologies, Inc.
203 Airport Road West
Belfast, Northern Ireland BT3 9ED
Tel: +44 28 90 93 94 00
Fax: +44 28 90 93 94 01
info@heartsine.co.uk



The products described in this brochure all meet the applicable European Medical Directive requirements.

CAUTION: U.S. Federal law restricts this device to sale by or on the order of a licensed practitioner.

© 2014 HeartSine Technologies, Inc. All rights reserved.
H009-014-122-0

Lifesaving, Pure and Simple

www.heartsine.com