

# HeartSine® Mobile AED Rescue Backpack

For samaritan® PAD Automated External Defibrillators



## Be Prepared for Cardiac Emergencies Even When You're on the Move

The HeartSine Mobile AED Rescue Backpack, custom made for the HeartSine samaritan® Public Access Defibrillator (PAD), lets you create and transport a heart safe solution to virtually any location.

Sudden cardiac arrest (SCA) can't always be prevented, but immediate access to an AED can increase the chance of survival by up to 74%, especially in remote areas where access to emergency services is often limited. In the event of an SCA, there's little time to react and even less time to think as an SCA can happen to anyone, at any time, regardless of age, gender or fitness.

The HeartSine Mobile AED Rescue Backpack has room for everything you will need to keep members of your group heart safe. This customized backpack can readily carry a HeartSine samaritan PAD AED with Pad-Pak™, a Paediatric-Pak™ or extra Pad-Pak, CPR prep kit and first aid kit.

### Ideal for Active Organisations

- Mountain/Water/Disaster Rescue Teams
- Community First Responders
- Scouting/Youth Groups
- Cycling, Hiking and Running Clubs
- School Athletics and Outings
- Family Outings

And more...



### Backpack Features

- Custom compartment for HeartSine AED with internal pocket for storage of Pad-Pak and Paediatric-Pak
- AED symbol for easy location of AED compartment
- First aid kit compartment with external symbol and small storage pocket
- D-rings for external attachment
- Large mesh pocket with adjustable cordlock to tighten
- Mesh side pockets with adjustable cordlock for storage of water bottle, cell phone and snacks
- Mesh shoulder strap pockets for flashlight, whistle or other small objects
- Adjustable shockcord with cordlock for keeping jacket within reach
- Adjustable chest and waist straps to fit adults or children
- Adjustable side straps to secure contents
- Interior mesh pocket to stabilize additional contents
- Padded back for extra comfort
- Adjustable padded shoulder straps
- Carrying handle
- Reflective safety strips



\* Backpack sold separately. HeartSine offers the AED, Pad-Pak, Paediatric-Pak and CPR Prep Kit.

## Why HeartSine AEDs?

- ▶ **Significantly smaller size** so it can travel with the group
- ▶ **Lightest AED** makes for easy transport in a backpack
- ▶ **Highest level of protection from dust and water (IP56)** offers unmatched robustness in both indoor and outdoor environments
- ▶ **Easy to use by all rescuers** due to simple one or two-button operation, intuitive voice/visual prompts and no lid to open nor complex displays or controls

## Backpack Specifications

- ▶ **Material:**  
Breathable polyester
- ▶ **Dimensions:**  
220 mm w x 360 mm h x 160 mm d  
(8.66 in w x 14.17 in h x 6.30 in d)
- ▶ **Weight:**  
Backpack only  
0.96 kg (1.9 lbs)  
  
Backpack with AED  
2.06 kg (4 lbs)
- ▶ **Colour:**  
Red with black trim



1. Valenzuela TD, et al. 2000. Outcomes of Rapid Defibrillation by Security Officers After Cardiac Arrest in Casinos. *New England Journal of Medicine*. 343:1206-09.

#### EMEA/APAC

HeartSine Technologies, Ltd.  
203 Airport Road West  
Belfast, Northern Ireland  
BT3 9ED  
Tel: +44 28 9093 9400  
Fax: +44 28 9093 9401  
info@heartsine.com



#### U.S./Americas

HeartSine Technologies LLC  
121 Friends Lane, Suite 400  
Newtown, PA 18940  
Toll Free: (866) 478 7463  
Tel: +1 215 860 8100  
Fax: +1 215 860 8192  
info@heartsine.com



The HeartSine products described in this brochure meet the European Medical Directive requirement.

UL Classified.  
See complete marking on product.

H009-032-337-3

The SAM 500P is not available for sale in the U.S.

© 2017 HeartSine Technologies LLC. All rights reserved.

[www.heartsine.com](http://www.heartsine.com)



HeartSine®