

HeartSine® samaritan® PAD Family



	HeartSine SAM 350P	HeartSine Connected SAM 350P	HeartSine SAM 360P	HeartSine Connected SAM 360P	HeartSine SAM 500P with CPR Advisor™	HeartSine Connected SAM 500P
ICG-based real-time feedback on CPR rate and force ("Push harder", "Push faster", "Push slower", "Good compressions")					X	X
CPR coaching with metronome	X	X	X	X	X	X
Integrated Wi-Fi connectivity for simplified readiness and remote AED program management		X		X		X
Communicates with LIFELINKcentral™ AED Program Manager to enable AEDs to be managed across a single or multiple locations		X		X		X
Fully automatic shock delivery			X	X		
Semi-automatic shock delivery	X	X			X	X
Highest level of protection from dust and water ingress in the industry (IP56)	X	X	X	X	X	X
Integrated battery and electrode for one-change maintenance	X	X	X	X	X	X
4-year electrode and battery life	X	X	X	X	X	X
Low energy escalating protocol (SCOPE™)	X	X	X	X	X	X
Operating impedance range of 20 – 230 ohms	X	X	X	X	X	X
Audible and visual metronome for CPR rate guidance	X	X	X	X	X	X
Compact, portable, and lightweight	X 1.1 kg (2.4 lb)	X 1.285 kg (2.83 lb)	X 1.1 kg (2.4 lb)	X 1.285 kg (2.83 lb)	X 1.1 kg (2.4 lb)	X 1.285 kg (2.83 lb)
Pediatric compatible*	X	X	X	X	X	X
90-minute ECG recording time	X	X	X	X	X	X
Automatic self test	X	X	X	X	X	X
Field upgradeable software	X	X	X	X	X	X
Carry case included	X	X	X	X	X	X

* For the SAM 500P, when used with Pediatric-Pak™, CPR Advisor is disabled.

For further information contact us at heartsinesupport@stryker.com or visit our website at heartsine.com.



HeartSine Technologies, Ltd.
203 Airport Road West
Belfast, BT3 9ED
United Kingdom
Tel: +44 28 9093 9400
Fax: +44 28 9093 9401



UL Classified. See complete marking on product.

The SAM 500P is not available for sale in the U.S.

© 2019 HeartSine Technologies. All rights reserved. H009-043-202-1 IE