


**HeartSine**<sup>®</sup>  
samaritan<sup>®</sup> PAD  
Family

Accessories

**stryker**



**Adult Pad-Pak**  
For patients > 8 years or  
25 kg (55 lb).

**Pad-Pak-03**  
**11516-000003**



**Pediatric-Pak™**  
For patients 1 to 8 years or  
up to 25 kg (55 lb).

**Pad-Pak-04**  
**11516-000004**



**Adult Pad-Pak**  
With TSO-C142a certification  
for use on aircraft.

**Pad-Pak-07**  
**11516-000027**



**USB Data Cable**  
For use with the Saver EVO™  
data management software  
to download event data and to  
perform software upgrades to  
the HeartSine samaritan PAD  
non-connected devices.

**PAD-ACC-02**  
**11516-000018**

## Cabinets and Mounting Options



**Rotaid Plus Wall Cabinet  
with Alarm**  
For indoor use.

**11996-000442** (Green)  
**11996-000441** (White)  
**11996-000443** (Red)  
**11996-000444** (Yellow)




**Rotaid Solid Plus Wall  
Cabinet with Alarm**  
For indoor or outdoor use.

**11996-000446** (Green)  
**11996-000445** (White)  
**11996-000447** (Red)  
**11996-000448** (Yellow)



**Rotaid Solid Plus Heat  
Wall Cabinet with Heat  
and Alarm**  
For outdoor use.

**11996-000450** (Green)  
**11996-000449** (White)  
**11996-000451** (Red)  
**11996-000452** (Yellow)



**AED Wall Bracket**  
(Metal)

Dimensions:  
15.6 H x 17.5 W x 11.4 D cm  
(6.14 H x 6.88 W x 4.48 D in)

**PAD-CAB-02**  
**11516-000023**

## Training Tools



**HeartSine samaritan PAD  
Trainer SAM 350P**

**TRN-350-XX\***  
**11516-000059**



**HeartSine samaritan PAD  
Trainer SAM 360P**

**TRN-360-XX\***  
**11516-000091**



**HeartSine samaritan PAD  
Trainer SAM 500P**

**TRN-500-XX\***  
**11516-000092**



**Trainer-Pak**  
(Replacement Electrode  
Cartridge)

**TRN-PAK-04**  
**11516-000017**

\* XX represents two-character country code.

# HeartSine Gateway™



**HeartSine Gateway**  
Includes removal tool, batteries  
and HeartSine Connected AED  
carrying case.  
**ACC-GTW-XX\***



**HeartSine Gateway  
Removal Tool**  
**ACC-REM-GW**



**HeartSine Gateway  
Batteries**  
**ACC-BAT-GW**

## Storage Options



**Wall Cabinet with Alarm,  
Custom Made**  
Dimensions:  
30.48 H x 27.94 W x 12.7 D cm  
(12 H x 11 W x 5 D in)  
Weight: 2.27 kg (5 lb)  
**PAD-CAB-04**  
**11516-000024**



**HeartSine AED Carrying  
Case**  
**PAD-BAG-01**  
**11516-000022**



**Mobile AED  
Rescue Backpack**  
**PAD-BAG-02**  
**11516-000114**



**Replacement electrodes**  
(for use with Trainer-Pak)  
**TRN-ACC-02** (10/pak)  
**11516-000009**  
**TRN-ACC-03** (25/pak)  
**11516-000011**



**HeartSine samaritan PAD  
Trainer Remote**  
**TRN-ACC-16** (SAM 350P)  
**11516-000014**  
**TRN-ACC-18** (SAM 360P)  
**11516-000016**




**TRN-ACC-17** (SAM 500P)  
**11516-000001**



**Battery Charger**  
(Replacement)  
**TRN-ACC-14**  
**11516-000012**

For further information contact us at [heartsinesupport@stryker.com](mailto:heartsinesupport@stryker.com) or visit our website at [heartsine.com](http://heartsine.com).

**EMEA/APAC**

HeartSine Technologies, Ltd.  
203 Airport Road West  
Belfast, Northern Ireland  
BT3 9ED  
Tel: +44 28 9093 9400  
Fax: +44 28 9093 9401 

**U.S./Americas**

HeartSine Technologies LLC  
121 Friends Lane, Suite 400  
Newtown, PA 18940  
Toll Free: (866) 478 7463  
Tel: +1 215 860 8100  
Fax: +1 215 860 8192

The HeartSine products described in this brochure meet the European Medical Directive requirement.



UL Classified. See complete marking on product.



H009-016-042-3 IE

© 2018 HeartSine Technologies LLC.  
All rights reserved.