

HeartSine® samaritan® PAD

Family



	HeartSine SAM 350P	HeartSine Connected SAM 350P	HeartSine SAM 360P	HeartSine Connected SAM 360P	HeartSine SAM 500P with CPR Advisor™	HeartSine Connected SAM 500P with CPR Advisor™
ICG-based real-time feedback on CPR rate and force ("Push harder", "Push faster", "Push slower", "Good compressions")					X	X
CPR coaching with metronome	X	X	X	X	X	X
Integrated Wi-Fi connectivity for simplified readiness and remote AED program management		X		X		X
Communicates with LIFELINKcentral AED Program Manager to enable AEDs to be managed across a single or multiple locations		X		X		X
Fully automatic shock delivery			X	X		
Semi-automatic shock delivery	X	X			X	X
High level of protection from dust and water ingress (IP56)	X	X	X	X	X	X
Integrated battery and electrode for one-change maintenance	X	X	X	X	X	X
4-year electrode and battery life	X	X	X	X	X	X
Low energy escalating protocol (SCOPE)	X	X	X	X	X	X
Operating impedance range of 20 – 230 ohms	X	X	X	X	X	X
Audible and visual metronome for CPR rate guidance	X	X	X	X	X	X
Compact, portable, and lightweight	X 1.1 kg (2.4 lb)	X 1.285 kg (2.83 lb)	X 1.1 kg (2.4 lb)	X 1.285 kg (2.83 lb)	X 1.1 kg (2.4 lb)	X 1.285 kg (2.83 lb)
Pediatric compatible	X	X	X	X	X*	X*
90-minute ECG recording time	X	X	X	X	X	X
Automatic self test	X	X	X	X	X	X
Field upgradeable software	X	X	X	X	X	X
Carry case included	X	X	X	X	X	X

* When used with Pediatric-Pak, CPR Advisor is disabled.

If you purchased your HeartSine Connected AED from an authorized Stryker distributor or reseller, this distributor or reseller will have access to your LIFELINKcentral AED Program Manager account and may receive notifications prompted by the HeartSine Connected AED. Please note that this setting to notify your distributor or reseller can be disabled at ANY time: if you wish to disable this setting, please send a request to Stryker Customer Support to self-manage your site without notifications to your distributor or reseller.

All claims valid as June 2020.

For further information, please contact your Stryker representative or visit our web site at strykeremergencycare.com

Emergency Care Public Access

AED users should be trained in CPR and in the use of the AED.

Although not everyone can be saved, studies show that early defibrillation can dramatically improve survival rates. AEDs are indicated for use on adults and children. AEDs may be used on children weighing less than 25 kg (55 lb) but some models require separate defibrillation electrodes.

The information presented is intended to demonstrate Stryker's product offerings. Refer to operating instructions for complete directions for use indications, contraindications, warnings, cautions, and potential adverse events, before using any of Stryker's products. Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your representative if you have questions about the availability of Stryker's products in your area. Specifications subject to change without notice.

Stryker or its affiliated entities own, use, or have applied for the following trademarks or service marks: CPR Advisor, HeartSine, LIFELINKcentral, Pediatric-Pak, samaritan, SCOPE, Stryker. All other trademarks are trademarks of their respective owners or holders. All other trademarks are trademarks of their respective owners or holders.

The absence of a product, feature, or service name, or logo from this list does not constitute a waiver of Stryker's trademark or other intellectual property rights concerning that name or logo.

06/2020

HeartSine SAM 500P is not available for sale in the U.S.

HeartSine Gateway is not available in all countries.

Made in U.K.

Date of Issue: 2020/06

H009-043-202-2

Copyright © 2020 HeartSine Technologies. All rights reserved.



HeartSine Technologies Ltd.
203 Airport Road West
Belfast, BT3 9ED
United Kingdom
Tel +44 28 9093 9400
Fax +44 28 9093 9401
heartlinesupport@stryker.com
heartline.com

Stryker European Operations B.V.
Herikerbergweg 110
1101 CM Amsterdam
Netherlands
Tel +31 (0)43 3620008
Fax +31 (0)43 3632001



UL Classified. See complete marking on product.