

HeartSine samaritan® PAD Family



	SAM 350P Semi-Automatic	SAM 360P Fully automatic	SAM 500P with CPR Advisor™
ICG-based real-time feedback on CPR rate and force ("Push harder", "Push faster", "Push slower", "Good compressions")			✓
Fully automatic shock delivery		✓	
Semi-automatic shock delivery	✓		✓
Highest level of protection from dust and water ingress in the industry (IP56)	✓	✓	✓
Integrated battery and electrode for one-change maintenance	✓	✓	✓
4-year electrode and battery life	✓	✓	✓
Low energy escalating protocol (SCOPE™)	✓	✓	✓
Operating impedance range of 20 – 230 ohms	✓	✓	✓
Audible and visual metronome for CPR rate guidance	✓	✓	✓
Compact, portable, and lightweight 1.1 kg (2.4 lbs)	✓	✓	✓
Paediatric compatible <i>(When used with Paediatric-Pak™, CPR Advisor is disabled)</i>	✓	✓	✓
90-minute ECG recording time	✓	✓	✓
Automatic self test	✓	✓	✓
Field upgradeable software	✓	✓	✓
Carry case included	✓	✓	✓

EMEA/APAC

HeartSine Technologies, Ltd.
203 Airport Road West
Belfast, Northern Ireland
BT3 9ED
Tel: +44 28 9093 9400
Fax: +44 28 9093 9401
info@heartsine.com



U.S./Americas

HeartSine Technologies LLC
121 Friends Lane, Suite 400
Newtown, PA 18940
Toll Free: (866) 478 7463
Tel: +1 215 860 8100
Fax: +1 215 860 8192
info@heartsine.com



The HeartSine products described in this brochure meet the European Medical Directive requirement.



UL Classified.
See complete marking on product.

H009-036-202-3 Intl

The SAM 500P is not available for sale in the U.S.
© 2017 HeartSine Technologies LLC. All rights reserved.

www.heartsine.com



HeartSine®